Airmatic™ II

- D** Gebrauchsanweisung
- UK** Operating instructions
- F** Mode d'emploi

Inhalt

Seite

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Table of contents

Page

1.0	Introduction	
1.1	Notes on operating instructions.....	5
1.2	Intended use	7
1.3	Function	7
1.4	Explanation of symbols	9
2.0	Safety advice	11
3.0	Installation and start-up	
3.1	Standard delivery	13
3.2	Illustrations	14 - 15
3.3	Connections	17
3.3.1	Electrical connection	17
3.3.2	Connecting a nystagmograph	17
3.3.3	Connecting a potential balance	17
3.3.4	Connection of the handle	17
3.3.5	Air inlet	17
3.4	Start-up	17
4.0	Operation	
4.1	Adjusting temperatures	19
4.2	Selecting temperature levels	19
4.3	Adjusting stimulation time	19
4.4	Description of operating modes	21
4.4.1	Preparation mode	21
4.4.2	Stimulation mode	23
4.4.3	Standby mode	23
5.0	Cleaning	
5.1	General information on cleaning and disinfection	25
5.2	Recommended instrument disinfectants	27
5.3	Recommended surface disinfectants	27
6.0	Maintenance and service	
6.1	Functional test	31
6.2	Visual inspection	32
6.3	Electrical test	32
6.4	Replacing of filter	32
7.0	Trouble shooting	36
8.0	Accessories and spare parts	
8.1	Accessories	39
8.2	Spare parts	39
9.0	Specifications	42
10.0	Disposal	44

General Standard Terms and Conditions

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1.0 Introduction

1.1 Notes on operating instructions

- The operating instructions contain important information on how to operate the Airmatic™ II safely, correctly and effectively. The operating instructions are intended for new operating personnel and also should be continued to be used as a reference manual. Reading and understanding the operating instructions will assist the user in avoiding risks, reducing repair costs, lowering down-time, increasing reliability and extending the service-life of the equipment. **The operating instructions must always be kept available and near the equipment.** Prior to first use, please read the chapter "Safety advice", to be prepared for any possible dangerous situations.

The basic principles are:

Judicious and careful work provides the best protection against accidents!

Operational safety and readiness for use depend not only on your capabilities, but also on care and maintenance given to the Airmatic™ II. Regular cleaning and service work are a must. Major maintenance and repair work should be carried out only by expert personnel authorised by the producer. For repairs, only use GN Otometrics spare parts to preserve the equipment warranty and ensure operational safety, readiness for work and to maintain the value of the equipment.

- The product Airmatic™ II bears CE Marking CE-0124 according to the EEC guideline of the council for medical products 93/42/EEC and meets the basic requirements of this guideline.
- The quality management system applied at the Producer has been certified according to international standards EN ISO 9001 and EN ISO 13485.
- Reprints, and text extracts are allowed only with written permission from GN Otometrics.

Abbreviations / symbols in these operating instructions:

- Indicating a list
 - Subdivision of a list/activity

The recommended sequence must be followed in each case!

Indicating particularly important advice!

↪ Describing the effect of an activity.

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UK

1.2 Intended use

- Compact unit for stimulation of the vestibule. Distinguished by its easy operation, it can be used in hospital and in a practice.
- This unit may not be used if the auditory canal is inflamed or contaminated.
- This unit is designed for short-time use on patients.

1.3 Function

- After activating the main switch, the optical readout displays tested.
- The unit automatically switches into the standby mode in which the heating and the pump are switched off.
- The unit can be changed into the stimulation mode for stimulating the vestibule. The Airmatic TM II is equipped with a timer for preselecting the stimulation time.

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UK

1.4 Explanation of symbols

Pay attention to operating instructions !
as to ISO /7000/0434 DIN 30600/1008
IEC 348

Type B equipment as to IEC 417

Fuse
as to IEC 417/5016, DIN 30600/0186

°C

Temperature in degree centigrade

S

Timer adjustment in seconds

Start

Stop

Timer

Cold stimulation level

Warm stimulation level

Heating and air flow ON

Heating and air flow OFF (standby)

Control output for connecting a nystagmograph
(graphical recorder as to DIN 30600,
IEC 417 5192)

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
⇒ Page 12

2.0 Safety advice

- The Airmatic™ II is produced according to IEC 601/EN 60601 and listed in the following classes:
 - VDE Class of protection 1
 - Class IIa (EEC 93/42).
- The equipment must only be connected to a power outlet with a correctly installed ground/earthing contact according to local and state regulations.
- Regulations for users in the European Economic Area: The Airmatic™ II must be subjected annually to a safety-technical inspection according to the regulations in these operating instructions.
- The Airmatic™ II may be used in supervised operation by qualified personnel that have been authorised by GN Otometrics and been trained for operating the equipment (IEC 601-1/EN 60601-1).
- The mains voltage specified on the type plate must match the power supply system.
- Prior to every use of the equipment, verify that the device is technically safe and in proper working condition. **Damaged leads** must be replaced immediately!
- Correct configuration in assembly of country-specific connections

green/yellow:	protective conductor (PE)
blue:	neutral conductor (N)
black resp. brown:	phase (L)
- The control panel must be visible and in reach of the operator. Ensure the equipment is stable prior to use.
- Prior to application, the air temperature must be checked by the user!
- Switch off the main power switch after using the equipment.
- The Airmatic™ II should only be operated in rooms used for medical purposes, and never in areas subject to explosion hazards.
- All additional equipment, which is connected to the analog and digital interfaces of the unit, must meet the requirements of relevant EN specifications (for inst. EN 60950 for data processing equipment and EN 60601 for electrical medical appliances). In addition, configurations must satisfy system specification EN 60601-1-1. When additional equipment is connected to the signal input or signal output section on the unit, the person making the connection is deemed "a system configuration operator" and as such is responsible for meeting the requirements of system specification EN 60601-1-1. For answers to additional questions, please contact your local specialist supplier or GN Otometrics technical service.
- GN Otometrics is not liable for personal injury and damage to property if
 - original GN Otometrics parts are not used,
 - the operating instructions are not completely followed,
 - assembly, new settings, alterations, extensions and repair have been carried out by personnel not authorised by GN Otometrics.
- The validity for the enclosed certificate of conformity expires if the customer or a third party manipulates the unit, e.g. modifications of any kind, by installing non-authorized accessories, removing warning or information labels as well as using the unit for inappropriate applications.

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UK

3.0 Installation and start-up**3.1 Standard delivery***Quantity*

1	Airmatic™ II Basic unit
1	Power cable
1	Handle with tube
1	Hose tips (30 pcs.)
1	Operating instructions

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D

Fig. 1. Front view

- ① Main switch
② Key switch f. heating/air flow ON/OFF (standby)
③ Coding switch for warm-stimulation level
④ Coding switch for cold-stimulation level
⑤ Coding switch for stimulation time
⑥ Key switch for selecting the warm-stimulation level (e.g. 44°C)
⑦ Key switch for selecting the cold-stimulation level (e.g. 30°C)
⑧ Key switch for start/stop of the stimulation
⑨ Temperature display (two-figure number, resolution 1°C), current value indication
⑩ Display of stimulation time (two-figure number, resolution 1s)
⑪ Support for handle
⑫ Connection for handle

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D

Fig. 2. Rear view

- ❶ Control output for a nystagmograph
- ❷ Unit plug with fuse compartment
- ❸ Equipotentiality connection
- ❹ Air inlet filter

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D

Fig. 3. Handle

- ❶ LED for indicating the stimulating process
- ❷ Timer Start/Stop key

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3.3 Connections

3.3.1 Electrical connection

- Connect the power cable to the socket on the unit (❷, fig. 2).
- Insert the power plug in a correctly installed socket with earthing contact.

3.3.2 Connecting a nystagmograph

- When controlling an ENG (electro-nystagmograph) or a CNG (computer-nystagmograph) at output (❶, fig. 2), please connect only recording equipment approved by GN Otometrics. Connecting cable available from GN Otometrics (see section 8.0).

3.3.3 Equipotentiality connection

- Connection for potential compensation (❸, fig. 2). Connecting cable available from GN Otometrics (see section 8.0).

3.3.4 Connection of the handle

- Only the handle which is intended for this purpose must be used:
 - Press the special connector plug slightly onto the socket in the unit.
 - Fasten it to the housing of the unit by turning the holding screws to the right.

Do not buckle the air tube !

3.3.5 Air inlet

- The ambient air is sucked into the unit through the filter.
- Screw connector for the air filter.

3.4 Start-up

- Insert the handle in its holder; the air outlet must point to reverse side of the unit.
- Switch on the unit (❶, fig. 1).
- Automatic display test with digital numbers "8 8" and audible signal.
- Automatic change into the standby mode.

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UK

4.0 Operation

- Use the main switch to start the unit.
- The unit is ready for operation after 1 second.

4.1 Adjusting temperatures

- Two variable temperature levels (20°C - 47°C) (48°C + 49°C only for testing purposes).

The lowest achievable stimulation temperature is approx. 2°C over the ambient temperature.

- Temperate setting by coding switches (③, ④, fig. 1).
 - left switch: for adjusting the "ten" partition
 - right switch: for adjusting the "one" partition
 - ↳ lower keys (+): Temperature increase
 - ↳ upper keys (-): Temperature decrease
- Standard settings:
 - Level for cold stimulation: 30°C
 - Level for warm stimulation: 44°C

4.2 Selecting temperature levels

- For selecting the desired temperature level, use the respective key (⑥, ⑦, fig. 1).
 - ↳ Display of the active level by illuminated diodes.
 - ↳ Display of the air temperature (current value) in °C.
- For switching off the heating system press the respective key (⑥, ⑦, fig. 1) of the active temperature level.
 - ↳ Light of LED of temperature level disappears.
 - ↳ Air temperature is shown (actual value) in °C.

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4.4.1 Preparation mode

Purpose:

The temperature set by the user is adjustable.

Properties:

- Temperature: corresponds to the preselected cold or warm stimulation level
- Air flow: 5.0 l/min.

Activation:

- By operating one of the key switches for selecting a temperature level (⑥, ⑦, fig. 1) or by operating the key switch for heating/air flow ON/OFF (②, fig. 1) in the standby mode
- If the active temperature push button is repeatedly pressed, the heating is switched off.
 - ↳ Air with a temperature that almost corresponds to the ambient air temperature is available.

Deactivation:

- Operating the key switch heating/air flow ON/OFF (②, fig. 1) switches the unit to the standby mode.
- Automatically switches to the standby mode when the unit is not used within a 3 minute period.

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Stimulation of the vestibule.

Properties:

- Temperature: Set according to pre-set cold or warm stimulation level
- Air flow: 5 l/min
- Duration: As pre-set by the timer.

Activation:

- Initially select the type of stimulation by actuating either warm-stimulation or cold-stimulation key (⑥, ⑦, fig. 1) (see section 4.1 for pre-setting of temperature).
- Operate the "timer-start-key" on the unit or on the handle.
- Preparation for stimulation:
 - As long as the push button is pressed, the pump remains switched off to allow the jet connection to be positioned in the auditory canal.
- When the push button is released, the thermal stimulation is performed for the time set by the user (⑤, fig. 1).
 - ↳ LED on handle (①, fig. 3) lights up during thermo-stimulation.
- At the end of the stimulation period a control signal for a recording unit is issued at the nystagmograph output.
- After completion of stimulation the pump is switched off.
- Repeated actuation of the "timer-start-key" during the stimulation leads to stopping timer operation.
- Second actuation of the currently active key effects deactivation of the pertaining stage.
 - ↳ Heating is switched off.
 - ↳ Stimulation with near-ambient air temperature.

4.4.3 Standby mode

Purpose:

- Reduction of energy consumption.
- Reduction of noise level.

Activation:

- Actuation of key "heating/pump on/off" (②, fig. 1).
 - ↳ Heating is switched off.
 - ↳ Pump is switched off after 2 s.
- Automatically switches to the standby mode after each stimulation process.
- Automatically switches to the standby mode when the unit is not used within a 3 minute period.

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UK

5.0 Cleaning

5.1 General information on cleaning and disinfection

After use, all parts which come into contact with the patient (jet connection, hose tips) must be removed and disinfected resp. autoclaved! You may use all cleansers and disinfectants listed in section 5.2 can be used.

The surfaces of the Airmatic™ II resist most of the common surface disinfectants.

However, do not use

- disinfectants containing concentrated organic or anorganic acids or bases, since these may cause corrosion damages;
- disinfectants containing chloramides, phenol derivatives or anionic tensides, since these may cause stress cracks in the plastic material used for the housing of the unit.

For cleaning and disinfection, you may also use disinfectant sprays or disinfectant tissues.

Switch off the unit using the main power switch before starting to clean and disinfect.

Wipe the unit surface with a cloth moistened with a cleaning or disinfecting solution. Make sure that no liquid penetrates the unit. You may use all cleansers and disinfectants listed in section 5.3.

Always observe the instructions for use by the manufacturer of the disinfectants, including all concentration specifications.

The above measures for cleaning and sterilizing do not replace the precautions for operating the unit.

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6.0 Maintenance and service

Change of filter (see section 6.4) should be carried out immediately by the user when required.

Safety-technical inspection

6.1 Functional test

Test

- ☐ ● Switch unit several times on and off.
 - ↳ Correct function of mains pilot lamp.
 - ↳ Automatic activation of displays test.
- ☐ ● Testing the displays.
 - ↳ All optical and audible indicators issue their signal during the first second after switching on (in particular segmental test "8 8 ").
- ☐ ● Set all temperature setting switches on 40°C and test consecutively by touching the corresponding keys whether in all ranges 40°C is reached in ACTUAL temperature indication.
 - Check air temperature.
 - Carry out the same test with temperatures 37°C and 28°C.
 - ↳ The particular combination of temperatures 40°C, 37°C and 28°C, as well as the test in all ranges ensures a comprehensive test of all contacts of the pre-setting switches, in ON- as well as in OFF-position, and also testing of all components involved in temperature control of the DESIRED temperature control system. Therefore, carrying out a test in one range only, or at only one temperature is not permissible.
- ☐ ● Set temperature at 49°C and close the air exhaust opening on the handle.
 - ↳ The electronic excess-temperature protection must respond after approx. 100-120 sec. (in the temperature display error code "F7" appears).
- ☐ ● Measure temperature of air at standard settings: 30° and 44°C.
 - ↳ The value may not deviate more than $\pm 0,5^\circ \pm 1/2$ digit, i.e. 1°C in total. Typically, however, values are better than $\pm 0,5^\circ$ C.

This measurement is critical !

Measuring probes or thermometers which are too large may cause faulty measurements because of substantial heat removal, inertia, and because their temperature sensitivity is not localised on a point but spread over a large surface. Use measuring instruments with an immersion probe, having a minimum precision of 0.3°C and a resolution of 0.1°C.

After setting the temperature, wait at least 1 minute before measuring to allow temperature stabilization of the double-barrelled hose and handle.

After setting the temperature 1 minute at least should be waited in order that the probe attains correct temperature.

6.2 Visual inspection

Test

- ☐ ● Check the unit for external damage (damage to cables, hoses, filter and handle).
- ☐ ● Check filter for dirt or a defect.
- ☐ ● Inspect unit, disconnected from the main power supply and with open housing, for internal damage, and leaks. Close housing again and tighten screws. Re-connect the main power cable.
- ☐ ● Check fuses for correct protection value and melting characteristics.
- ☐ ● Check sign plates and inscriptions for completeness and readability.

6.3 Electrical test

Test

- ☐ ● Measure protective conductor resistance according to VDE 0751 (max. 0.3 Ohm).
- ☐ ● Measure and compare with "first value measured" of unit leakage current according to VDE 0751 (measured during operation).
 - ↳ The "value measured first" may not be exceeded by more than 50%.
- ☐ ● Document and record the evaluation and tests according to VDE 0751.

6.4 Replacing the filter

- The air filter (④, fig. 2) must be replaced once per year as follows:
 - Remove the filter by unscrewing it from the holder (turn the filter to the left).
 - Mount the new filter by screwing it into the holder (turn the filter to the right).

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

7.0 Trouble-shooting

If, errors cannot be corrected inform the GN Otometrics service staff. Do not attempt to repair the unit yourself!

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8.0 Accessories and spare parts

8.1 Accessories

<u>Description</u>	<u>Article-No.</u>
Air inlet filter	000.0545.0
Power cable	507.0859.0
Handle, complete	502.1035.0
Hose tips for jet (30 pcs.)	502.0844.0
Operating instructions	502.1100.i

8.2 Spare parts

<u>Description</u>	<u>Article-No.</u>
Unit basis, to be plugged in	000.0639.0
Fuse T800mA H / 250V 5x20mm	008.0081.0
Air inlet filter	000.0545.0
Power cable	507.0859.0
Handle, complete	502.1035.0
Jet connection	502.1045.0
Hose tips for jet (30 pcs.)	502.0844.0

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F

⇒ Page 43

Voltage range	100 - 240 V~ ± 5 %; 50/60 Hz
Current input	0.75 A max.
Power cable	3 m
Power consumption	85 W max.
Operating modes	Preparation mode; Stimulation mode (at temperature preselected for the cold stimulation level resp. warm stimulation level); Heating off and no air flow (economy mode, standby mode)
Stimulation time	can be preselected by timer from 1 up to 99 sec.
Timer indication	Accuracy ± 0.5 s ± ½ digit
Air temperature	20°C - 47°C
Lowest temperature	approx. 2°C above room temperature
Temperature indication	Accuracy ± 0.5 °C ± ½ digit
Temperature deviation	< ± 1°C
Air flow	5.0 l/min ± 10 %
Connections	Power connection; Control output for a nystagmograph; Equipotential equalization; Connection for the handle; Air inlet
Operating time	Short-time operation: 1. Automatic shut-off after completion of stimulation. 2. Automatic shut-off after 3 minutes.
Fuses	quick-acting 8 A (220-240 V~) quick-acting 15 A (110-127 V~)
Protective earth conductor resistance	0.1 Ω max.
Earth leakage current	0.5 mA max.
Enclosure leakage current	0.1 mA max.
Patient leakage current	0.1 mA max.
Ambient conditions	
Transport/storage	-20...+50°C; 5...90 % humidity without condensing, air pressure 700...1060 hPa
Operation	+10...+35°C; 20...80 % humidity without condensing, air pressure 700...1060 hPa
Dimensions HxWxD	14.5 x 37 x 32 cm
Weight	approx. 3.7 kg
Protection class (EN 60601-1)	I
Degree of protection	Type B 
Protection category	IPX0
Classification acc. to Annex IX EEC directions 93/42/EEC	IIa
CE marking	CE 0124
Rules applied	EN 60601-1: 1990 + A1:1993 + A2:1995 EN 60601-1-2: 1993 (EMV / EMC)
UMDNS-Code	10-548
REF	502.1100.0
Canadian Classification	
Device group	Ear, nose, throat
PNC	77ETP
Risk Class	2
Description	STIMULATOR, CALORIC (WATER)

Classification canadienne	
Application	ORL
PNC	77ETP
Classe de risque	2
Description	Stimulation calorique (eau)

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UK

10.0 Disposal

- Packaging material, cardboard and/or PE foam, can be fully recycled or returned to your supplier.
- The Airmatic™ II does not contain any hazardous materials.
- The material of the housing can be completely recycled.
- The component parts of the Airmatic™ II must be correctly discarded and the materials should be carefully separated for recycling.
- Electronics circuit boards must be discarded using the appropriate recycling process.
- Used hose tips, which no longer can be disinfected, must be discarded into domestic waste immediately.

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**EC - DECLARATION OF CONFORMITY
FOR MEDICAL PRODUCTS**

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GmbH & Co. KG